


**LABQUEST**  
BY **BOROSIL®**

# KJELDAHL DISTILLATION

**OPERATING MANUAL**

**KDI 030**

**KDI 040**



**DEALER :**



## THANK YOU NOTE

*We Borosil, one of India's most customer oriented brands truly appreciate your business and express our gratitude for the trust you have placed on us.*

*We hope your choice serves you well in your scientific endeavors and aspire to have the pleasure of doing business with you for years to come.*



**TABLE OF CONTENTS**

<b>Sr. No.</b>	<b>Particular</b>	<b>Page No.</b>
1.	Specifications of the product	6
2.	Packing list	7
3.	Safety and Warning	8
4.	Unboxing of the product	9
5.	Product Identification	10
6.	Dosing reagent bottle connection	15
7.	Installation of unit	17
8.	Tube Connection for the unit	18
9.	Operating the unit	20
10.	Errors	26
11.	Troubleshooting	27
12.	Warranty Registration	29
13.	Statement of Warranty	31
14.	Contact information	32

## SPECIFICATIONS

<b>PARAMETERS</b>	<b>KDI030</b>	<b>KDI040</b>
<b>Test tube volume</b>	250ml	250ml
<b>External dimensions (WxDxH)mm</b>	350x380x725	350x380x725
<b>Mains voltage</b>	220v-240v	220v-240v
<b>Current consumption</b>	4.8A	4.8A
<b>Power consumptions</b>	1000W	1000W
<b>Heater</b>	Coil heater	Coil heater
<b>Body material</b>	MS powder coated	MS powder coated
<b>Bottom tray material</b>	SS304	SS304
<b>Display</b>	LCD	LCD
<b>Temperature probe</b>	PT-100	PT-100
<b>Door close sensor</b>	Yes	Yes
<b>Test tube detection</b>	Yes	Yes
<b>Steam power</b>	25% to100%	25% to100%
<b>Dosing funnel</b>	NA	NA
<b>Glassware material</b>	Borosilicate glass	Borosilicate glass
<b>Max. Distillation time (In min/sample)</b>	2-9 min	2-9 min
<b>Frequency</b>	50-60 Hz	50-60 Hz
<b>Memory</b>	50 user programmable program	50 user programmable program
<b>USB port</b>	Yes	Yes
<b>Safety sensors</b>	Door sensor,Test tube sensor	Door sensor,Test tube sensor
<b>Sensor cooling water flow</b>	Yes	Yes
<b>Ambient temperature</b>	15°C to 40°C	15°C to 40°C
<b>Measuring range</b>	0.1mg to 200mg	0.1mg to 200mg
<b>Reproducibility</b>	<1%	<1%
<b>Recovery</b>	>99.5%	>99.5%

**PACKING LIST - BOX 1**

1.	Kjeldahl Distillation Unit.....	01 No.
2.	Power cable 5 Amps .....	01 No.
3.	Key for door .....	02 Nos.
4.	Leveling Rubber Shoe.....	04 Nos.
5.	Selection Knob .....	01 No.
6.	Drain Valve .....	01 No.
7.	T-Brass connector .....	01 No.
8.	S.S. Hose clamps.....	06 Nos.
9.	Silicone tube.....	10 mtr
10.	Braided hose pipe .....	04 mtr
11.	Operating Manual.....	01 No.
12.	Test Report.....	01 No.
13.	Cable tie 100mm .....	05 Nos.
14.	Drip Tray LH.....	01 No.
15.	Drip Tray RH.....	01 No.
16.	Brass connector male.....	01 No.
17.	Ammonium sulphate sample bottle 3ml .....	01 No.
18.	Oven mitts .....	01 No.

**PACKING LIST - BOX 2**

1.	Condenser .....	01 No.
2.	Dosing Head & Dosing holder.....	01 No. each
3.	Silicone Tubing for Dosing ID100mm length .....	02 Nos.
4.	Test tube (250 ml).....	02 Nos.
5.	Teflon Tube ID 8mm x 440mm .....	01 No.
6.	Silicon Tube for condenser 8mm x 210mm.....	01 No.
7.	Teflon Tube ID 8mm x 440mm .....	01 No.

**PACKING LIST - BOX 3**

1.	Aspiration bottle with level sensor	
	• KDI030.....	02 Nos.
	• KDI040.....	03 Nos.

## SAFETY AND WARNING



### Important operating and maintenance instructions

Read the accompanying text carefully.



### Potential electrical hazards

Only qualified persons should perform procedures associated with this symbol.



### Lifting hazard warning

- The Kjeldahl distillation unit weighs more than 30 kg. Take adequate safety measures when moving this device.
- Equipment being maintained or serviced must be turned off and locked off to prevent possible injury.



### CAUTION :

- Always use the proper protective equipments (clothing, gloves, goggles and face mask etc.)
- Always follow the good hygiene practices.
- Each individual is responsible for his or her own safety.
- Do not keep the unit on wet bench top.
- Do not wash the unit with water.
- Wipe the unit with dry cloth after usage.
- Ensure that no fluid is spilled on the top surface of the equipment.
- The unit should be plugged to standard 230V, 50HZ, 6A, 3 pin power socket.
- The unit should be plugged into the power socket, having proper earthing.
- Use only distilled water to operate the unit.



### WARNING :

Inadequate earthing at the installation facility can lead to hazardous electrical shocks. The manufacturer is not liable for any injury or death resulting from electrical hazards due to faulty earthing in the lab.

## UNBOXING OF THE PRODUCT

### UNPACKAGING INSTRUCTION

1. Remove the unit from the Box 1.

### FOR ACCESSORIES :

1. Remove from the pouch.


- Knob
- Door keys
- Hose pipe
- Power cable - 6 Amp
- Leveling Shoe
- Brass Tee Connector
- Drain Valve for Boiler

2. Remove from box 2.

- Condenser
- Dosing Head with teflon tube
- Silicone tubing 10 ID and 8 ID
- Test tube (300 ml)
- Teflon Tube

3. Remove from box 3.

- Remove Aspiration Bottles



## PRODUCT IDENTIFICATION - KDI030

- |                          |                              |
|--------------------------|------------------------------|
| A. DOSING HEAD           | H. CONDENSER WATER DRAIN     |
| B. DISTILLED WATER INLET | I. CONDENSER                 |
| C. NaOH INLET            | J. CONDENSER WATER INLET     |
| D. DOSING HEAD HOLDER    | K. CONDENSER COLLECTION TUBE |
| E. SAMPLE TEST TUBE      | L. CONICAL FLASK 250 ml      |
| F. STEAM INLET TUBE      | M. DRIP TRAY LH & RH         |
| G. TEST TUBE HOLDER      |                              |




Fig. No. 10.1

## PRODUCT IDENTIFICATION - KDI040

- |                            |                              |
|----------------------------|------------------------------|
| A. DOSING HEAD             | H. CONDENSER WATER DRAIN     |
| B. DISTILLED WATER INLET   | I. CONDENSER                 |
| C. NaOH INLET              | J. CONDENSER WATER INLET     |
| D. DOSING HEAD HOLDER      | K. CONDENSER COLLECTION TUBE |
| E. SAMPLE TEST TUBE        | L. BORIC ACID INLET TUBE     |
| F. STEAM INLET TUBE        | M. CONICAL FLASK 250 ml      |
| G. SAMPLE TEST TUBE LOADER | N. DRIP TRAY LH & RH         |




Fig. No. 11.1

## PRODUCT IDENTIFICATION

- A. DISPLAY PCB
- B. LCD DISPLAY
- C. SELECTION KNOB
- D. FRONT ACRYLIC DOOR
- E. DOOR HANDLE
- F. LEVELING SHOE
- G. B-TYPE USB SOCKET
- H. CONTROLLER PCB
- I. POWER PLUG SOCKET with FUSE (5 Amp)
- J. POWER SWITCH
- K. DISTILLED WATER INLET FOR BOILER
- L. STEAM OUT
- M. COOLING WATER INLET
- N. BOILER DRAIN
- O. EXCESS STEAM OUTLET
- P. COOLING WATER OUTLET




Fig. No. 12.1




Fig. No. 12.2




Fig. No. 12.3

## PRODUCT IDENTIFICATION - KDI030

- A. SOLENOID VALVE FOR STEAM IN
- B. DISTILLED WATER INLET FOR BOILER
- C. SOLENOID VALVE FOR BOILER FILLING
- D. BOILER WITH HEATER
- E. STEAM INLET FOR DOSING
- F. FLOW SWITCH (COOLING WATER)
- G. SOLENOID VALVE FOR (COOLING WATER)
- H. BOILER DRAIN VALVE




Fig. No. 13.1

## PRODUCT IDENTIFICATION - KDI040

- A. SOLENOID VALVE FOR STEAM IN
- B. DISTILLED WATER INLET FOR BOILER
- C. SOLENOID VALVE FOR BOILER FILLING
- D. BOILER WITH HEATER
- E. STEAM INLET FOR DOSING
- F. FLOW SWITCH (COOLING WATER)
- G. SOLENOID VALVE FOR (COOLING WATER)
- H. BOILER DRAIN VALVE




Fig. No. 13.2

## PRODUCT IDENTIFICATION




Fig. 14.1

Dosing Head

Product code :- BLG0DOH000KJLDIST




Fig. 14.2

Condenser

Product code :- BLG0CON000KJLDIST




Fig. 14.3

Dosing Head Holder

Product code :- BLM0DHH0000KD100V4




Fig. 14.4

Flow Switch

Product code :- BLMOFLS10000000KD1




Fig. 14.5

Test Tube Loader

Product code :- BLGAFUT010PTFEDI1



Fig. 14.6

Test Tube

Product code :- BLGATTQ01041250290

## DOSING REAGENT BOTTLE CONNECTION - KDI030




Fig. 15.1

1. Connect Bottle 1 level sensor wire to BOTTLE 1.
2. Connect Bottle 2 level sensor wire to BOTTLE 2.




Fig. 15.2

## DOSING REAGENT BOTTLE CONNECTION - KDI040





Fig. 16.1

(FOR SOIL TESTING)



1. Connect Bottle 1 level sensor wire to BOTTLE 1.
2. Connect Bottle 2 level sensor wire to BOTTLE 2.
3. Connect Bottle 3 level sensor wire to BOTTLE 3.




Fig. 16.2

## INSTALLATION OF UNIT

### CHECK BELOW POINTS BEFORE INSTALLATION OF THE UNIT.

1. The required water pressure for the unit should be around 5-6bar (70psi)
2. TDS value should be between 300ppm-500ppm, if the TDS value is more than 1200ppm, it is recommend to install the water softener or the appliance protection filter to/for the unit.
3. It is recommend to connect chiller to condenser unit for better results.
4. Distilled water is recommend for boiler to prevent scaling.
5. Boiler should be drained frequently after completing 10-20 cycles.

### INSTALLATION OF DOSING HEAD AND CONDENSER

1. Remove dosing head and condenser from box - 2.
2. Fix them with the unit as shown in the figure.
3. Connect Condenser water drain & Condenser water inlet with silicone tubing from the unit (Ref. fig. 13.1).
4. Connect Chemical dosing tubes (Ref. fig. 13.1) & Steam inlet tube to Dosing Head. (Ref. fig. 20.2)
5. Connect Boiler water inlet with hose pipe for operating the unit.
6. For connection of Dosing reagent bottles, insert tubings in their respective cans with level sensor connection on to the unit.  
(Ref page no. 17 &18 )




Fig. No. 17.1



## INSTALLATION WITH DOSING HEAD AND CONDENSER




Fig. 19.1

## OPERATING THE UNIT




fig. no.20.1

### • Main Menu (fig. no.20.1)

1. Rotate the knob clockwise and anti-clockwise for menu scroll.
2. Cursor on the left shows the position.
3. To select the particular menu "SINGLE PRESS" the knob.

SELECTION KNOB  
ROTATE KNOB TO SCROLL  
SINGLE PRESS TO SELECT  
LONG PRESS TO EXIT

### • START Menu (fig. no.20.2)

1. "SINGLE PRESS" to scroll down through the parameters, rotate the knob clockwise to increase the value of the parameter and rotate the knob counterclockwise to decrease the value of the parameter.
2. The arrow symbol ">" in the left of the display indicates the cursor position. "SINGLE PRESS" the knob after setting the last parameter value to START the distillation process.
3. "LONG PRESS" the knob to go to the first parameter of the window. "LONG PRESS" when on the first parameter to return back to the MAIN MENU window.




fig. no.20.2




fig. no.20.3

### • SETUP Menu (fig. no.20.3)

1. Rotate the knob clockwise and anti-clockwise for menu scroll.
2. Cursor on the left shows the position.
3. To select the particular menu "SINGLE PRESS".

• **PRESET MENU** (fig. no.21.1)




fig. no.21.1

1. There are 50 programmable programs from P0 to P49. The symbol “ > ” indicates the cursor. User can edit and save the programs as per their sample requirements.
2. “SINGLE PRESS” the knob to select the program to run or edit it. Rotate the knob clockwise and anticlockwise to scroll through the preset menu programs.
3. After selecting a particular program user can RUN or EDIT the program.


i) **RUN** -

User will be directed to the user input menu where all the parameters values will be taken as of the program selected. “SINGLE PRESS” the knob to start the distillation process. To edit the existing parameter values “LONG PRESS” the knob. User can see the cursor “>” on the first parameter and can change the values for the same. “SINGLE PRESS” & rotate the knob to scroll the parameters.


ii) **EDIT** -

User can edit the parameters value and save the same value in EDIT function. After selecting the edit option the parameter value which is seen on the display above will blink and user can identify that he/she is on which parameter. “SINGLE PRESS” the knob to scroll through the parameters and rotate the knob to change the values of the parameters. “SINGLE PRESS” the knob after setting the last parameter value i.e. SP (Steam Power) to save the changes made.

• **PUMP PRIMING** (fig. no.22.1)



1<sup>st</sup> STAGE



2<sup>nd</sup> STAGE

fig. no.22.1

1. Load test tube and conical flask before pump priming.
2. “SINGLE PRESS” to start pump priming.
3. User have to visually see the chemical coming in the test tube by “SINGLE PRESS” of the knob.
4. “SINGLE PRESS” every time will switch on the pump for priming.
5. “SINGLE PRESS” the knob to prime the next pump, then the previous pump will stop priming.

• **STEAM POWER** (fig. no.22.2)




fig. no.22.2

1. Rotate the knob to shift from one % power to another % (use cursor to locate the position)
2. To select a particular steam power between 25,50,75,100 rotate the knob & take the cursor below that value and “SINGLE PRESS”.
3. Once set, the user will not have to set it again unless the user wants to change the steam power.
4. “LONG PRESS” for auto return.

• **PUMPING TIME (FOR SERVICE ENGINEER)**

(fig. no.22.3)

1. This has to set to 2 sec.
2. One can change the value by rotating the knob.
3. To save this value “SINGLE PRESS”




fig. no.22.3

• **SENSOR STATUS** (fig. no.23.1)




fig. no.23.1

1. "SINGLE PRESS" the knob on SENSOR STATUS parameter in SETUP menu to select the setting. The arrow on the leftmost side of the screen indicates the cursor position. "SINGLE PRESS" to scroll down through the parameters. Here are two states in every parameter.

- i. **Override** - To override the sensor or to disable the sensor.
- ii. **Active** - To enable the sensor.

3. To save the changes done in sensor status "SINGLE PRESS" the knob when the cursor is on the last parameter, the changes made will be saved and will be retained even if the Instrument is powered OFF and ON.
4. To exit the setting without saving "LONG PRESS" the knob. It will exit and come to the previous screen / menu.

• **BOILER FILLING** (fig. no.23.2)

1. This indicates that the boiler is filling.
2. If there is no water flow at the inlet it will show the error of "INCREASE THE FLOW"




fig. no.23.2

• **HEATER ON** (fig. no.23.3)




fig. no.23.3

1. The heater is turned on.
2. User can see the current temperature of boiler on the screen.
3. The heater will remain ON until the temperature of the boiler reaches to 98°C .

• **DOSING REAGENT** (fig. no.24.1)

1. This window indicates the automatic dosing of acid, alkali and water for KDI040 .
2. It indicates automatic dosing of alkali and water for KDI030 .




fig. no.24.1

• **PROCESS** (fig. no.24.2)

1. This window shows the starting of steam dosing and distillation process.
2. It will display the input values provided by the user and dosing time.
3. The dosing time is the indication of remaining time for the process to be completed.
4. When flow rate is reduced in condenser error is displayed.(INCREASE FLOW RATE)




fig. no.24.2

**NOTE :**

- User have to increase the flow within 60 sec after display of message showing flow error.
  - If user fails to increase the flow, the process will be automatically terminated after 60 sec.
5. After completion of steam dosing it will return to the main menu on the display and user can remove the conical flask and the test tube.
  6. For next sample of distillation load the next test tube on loader and place the conical flask.
  7. After placing the test tube, single press the selection knob it will automatically start the next distillation process.

• **AFTER COMPLETION OF STEAM DOSING**

1. Remove the conical flask and "SINGLE PRESS" the knob to end the process. (fig. no.25.1)




fig. no.25.1




fig. no.25.2

2. This window indicates the completion of steam dosing process, with buzzer indication. (fig. no.25.2)

3. Remove the test tube when this window appears. (fig. no.25.3)




fig. no.25.3




fig. no.25.4

4. "SINGLE PRESS" the knob to continue next process.
5. "DOUBLE PRESS" the knob to start the cleaning cycle.
6. The cleaning cycle is of 4 min. (fig. no.25.4)

## ERRORS

- **DOOR CHECK** (fig. no.26.1)

1. This error window will appear when the door is open during the process.
2. This window will disappear after closing the door and the process will be resumed.




fig. no.26.1




fig. no.26.2

- **TUBE CHECK** (fig. no.26.2)

1. This error window will appear when test tube is not loaded into the loader before the process starts.
2. In process test tube removed.

- **FLOW CHECK** (fig. no.26.3)

1. This error window will appear when the flow rate is less in the condenser during the process & before the process starts.




fig. no.26.3




fig. no.26.4

- **LEVEL CHECK (ACID, ALKALI, WATER)**

(fig. no.26.4)

1. This error window will appear when the level of reagents is less in the reagent bottle.
2. This window will disappear after filling the bottle with reagent ( "SINGLE PRESS" the knob).

## TROUBLESHOOTING

### 1. If the unit is not turning ON.

- Check the power supply in AC mains.
- Make sure power cable is inserted to the socket properly.
- Check whether the main switch is ON or OFF.
- Check if the illuminated switch is OFF.
- Ensure the main switch is ON.

### 2. If the fuse is blown.

- Remove the power cable from unit.
- Remove the holder from the back side of the equipment, in the control panel box .
- Check the fuse, if it is damaged please change the fuse.





### WARRANTY REGISTRATION

Please handover this registration form to the distributor from where you have purchased this product. The warranty is valid only when this warranty registration form is received by us within 30 days from the date of purchase.

Product : KDI030  KDI040

Product Sr. No.: \_\_\_\_\_

Date of Invoice : \_\_\_\_\_

Invoice No.: \_\_\_\_\_

#### Customer name & address

Name : \_\_\_\_\_

\_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Telephone: \_\_\_\_\_

E-mail: \_\_\_\_\_

#### Customer sign & seal

#### Dealer name & address

Name : \_\_\_\_\_

\_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Telephone: \_\_\_\_\_

E-mail: \_\_\_\_\_

#### Dealer sign & seal



# **BOROSIL®** Scientific

STATEMENT OF WARRANTY

Borosil confirms that this product has been manufactured in accordance with our technical specifications and quality requirements.

- Borosil warrants the product from manufacturing and workmanship defects for a period of 12 months from the date of invoice.
- Warranty is void, if equipment and apparatus is not operated as prescribed in the operating manual supplied along with the unit.
- To be covered under warranty.
  - » Units have to be connected to standard 230V 50Hz, 6A wall sockets with proper earthing.
  - » Warranty does not cover rust and physical damage to metal parts due to corrosive environment in the lab and breakage of glass parts.

**Terms :**

- » In the event of malfunction due to defect, the buyer will have to follow the Borosil's service process.
- » Certain units can not be serviced / rectified at the buyer's place and the units may have to be brought to Borosil's service center as advised by Borosil's representatives.
- » In no event shall Borosil be liable for consequential or incidental damages.

<b>INVOICE DATE</b>	<b>BUYER</b>	<b>AFFIX SERIAL NUMBER</b>
<b>INVOICE#</b>		
<b>Dealer name &amp; address</b>		<b>Dealer sign &amp; seal</b>

**BOROSIL SCIENTIFIC LIMITED**

Corporate Office : 1101, Crescenzo G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai-400051, India



**: MANUFACTURED BY :**

**Borosil Scientific Limited**

Plot No. 7, Sr. No. 234, 235 & 245,  
Indialand Global Industrial Park,  
Hinjewadi Phase 1, Pune - 411057.

*Write to us on above address.*

**: MARKETED BY :**

**Borosil Scientific Limited**

1101, G-Block, Parinee Crescenzo,  
BKC, Bandra East, Mumbai - 51

Maharashtra, India.

**: CUSTOMER CARE CONTACT :**

**Phone : 1800 22 4551 | Email : [lab.support@borosil.com](mailto:lab.support@borosil.com)**

**Website : [www.borosilscientific.com](http://www.borosilscientific.com)**